

A.M.I. FINGER DISPOSABLE INSTRUMENT








Order Code EFI3001



A.M.I. EasyFinger Disposable
Disposable tip for A.M.I.
EasyInstrument Premium
Suitable for blunt dissection on the posterior side of organs and structures, and for positioning surgical loops or implants (e.g. for A.M.I. Soft Gastric Band, for A.M.I. B-Band, for fundoplication, etc.).
Sterile packed.
10 finger / box

A.M.I. AIGNER GASTRIC BAND CUTTER

Order Code	A.M.I. Aigner Gastric Band Cutter System	Technical Details
GBC 360 	The system consists of: _ GBC-Handle _ GBC-Toggle _ GBC-Toggle Holder _ GBC-Flexible Wire Guide _ GBC-Cutting Wires _ Using Manual	Delivered non sterile. Stainless steel. Fully detachable.
GBC 363 	A.M.I. GBC-Flexible Wire Guide Flexible wire guide to bring in the Cutting Wires. To be used in combination with the GBC 360. 1 Unit.	Stainless steel. Delivered non sterile.
GBC 364 	A.M.I. GBC-Cutting Wire For single use in combination with the GBC 360. Box of 3 Units.	Stainless steel. Delivered non sterile.
GBC 365 	A.M.I. GBC-Toggle Holder Safety spring for the GBC 360. 1 Unit.	Stainless steel. Delivered non sterile.



Gastric Band Cutter Migrasyona uğramış Gastrik Bantların cerrahi gerektirmeden endoskopik yöntemle çıkartılması için kullanılır.