

Mesh Sabitleme Sistemi



2D Reinforcement Prosthesis Microval Intra

MICROVAL INTRA **çift taraflı kompozit protez** olup, intra-peritoneal yaklaşımlarda, evantrasyon veya umbilical hernia tamirinde kullanılmaktadır. 70GSM densitivlik 100% nonwoven polpropylen, bir tarafı silikon ile kaplanmış, yapışıklık önleyici meshdir.

Intra abdominal hernia, open abdomen, diyagrafmatik hernia, insizyonel hernia ve batının kapatılamadığı durumlarda da kullanılmaktadır.

Paketleme: PET/PE - TYVEK çift steril, Karton kutu


Sterilizasyon: Ethylene Oxyd.

Référence	Désignation		
472012	MICROVAL-INTRA reinforcement prosthesis, Diameter 12 cm		
472015	MICROVAL-INTRA reinforcement prosthesis, rectangle 10x15 cm		
472517	MICROVAL-INTRA reinforcement prosthesis, rectangle 15x17 cm		
472522	MICROVAL-INTRA reinforcement prosthesis, cask shape, 10/15x22 cm		
472030	MICROVAL-INTRA reinforcement prosthesis, rectangle 20x30 cm		
472872	MICROVAL-INTRA reinforcement prosthesis, for hiatal hernia repair 7x8cm		

MICROVAL-3D-ANATOMIQUE PROSTHESIS

Reinforcement Prosthesis
for hernia repair in knitted
polypropylene

%100 biocompatible meshdir.



2D REINFORCEMENT PROSTHESIS IN KNITTED POLYPROPYLENE PRECUT SHAPE			
Référence	DESIGNATION		
413410	2D reinforcement prosthesis in knitted Polypropylene 10x14 cm for Laparoscopy surgery		
413612	2D reinforcement prosthesis in knitted Polypropylene 12x16cm for Laparoscopy surgery		
413105	2D reinforcement prosthesis in knitted Polypropylene 5x10cm for Lichenstein technicalsurgery		
413126	2D reinforcement prosthesis in knitted Polypropylene 6x12cm for Lichenstein technical surgery		
413127	2D reinforcement prosthesis in knitted Polypropylene 7x12cm for Lichenstein technical surgery		
413130	2D reinforcement prosthesis in knitted Polypropylene 9x13cm for Lichenstein technical surgery		

KNITTED POLYPROPYLENE MODUCLAR PLUGING SYSTEM			
Référence	DESCRIPTION		
411009	Modular plug 9cm diameter		
413004	Modular plug 9cm diameter + LF type Prosthesis with overlap 9x13 cm		
413006	Modular plug 9cm diameter + LF type Prosthesis with overlap 7 x 12 cm		
413008	Modular plug 9cm diameter + LF type Prosthesis with overlap 5 x 10 cm		
413009	Modular plug 9cm diameter + LF type Prosthesis with overlap 6 x 12 cm		
411126	PB type Mesh/Plug with prothesis 6x12cm and Plug 6 cm diameter		
413005	Plug disk diameter 5cm and type Lichtenstein prothesis 5x10 cm		
413007	Plug disk diameter 7cm and type Lichtenstein prothesis 9 x 13 cm		

ANATOMICAL AND PRE SHAPED : PROSTHESIS FOR LAPORSCOPY SURGERY			
Reference	Description		
411141	Pre-shaped prosthesis 10 x 14 cm		
411000	3D anatomical mesh small right 10x14 cm		
411001	3D anatomical mesh small left 10x14 cm		
411010	3D anatomical mesh medium right 12x16 cm		
411011	3D anatomical mesh medium left 12x16 cm		