

# OTSC SYSTEM FOR FLEXIBLE ENDOSCOPY

o o o o o o o o v e s c o



## OTSC<sup>®</sup> system....


Over-The-Scope-Clipping system for flexible endoscopy



# OTSC SYSTEM FOR FLEXIBLE ENDOSCOPY


## OTSC® TWIN GRASPER®

For approximating lesions before clip application



## OTSC® ANCHOR

For aligning the OTSC® applicator cap with the target tissue




### OTSC Anchor:

- Support treatment of hemorrhage
- NiTi anchor
- Anchors into tissue before clip release

### OTSC Anchor:

- Support treatment of hemorrhage
- NiTi anchor
- Anchors into tissue before clip release




...targeting the lesion...



...aligning OTSC applicator



...releasing the clip...



...result

## OTSC® CLIP



OTSC® traumatic clip

OTSC® atraumatic clip

### OTSC clip:

- Shape memory alloy
- Mounted to scope tip
- Unique tissue approximation & copression

## OTSC® HAND WHEEL



### OTSC system:

- Ancillary equipment
- Tissue manipulation devices
- Disposable



## OTSC® APPLICATOR CAP

## OTSC® ANCHOR

## OTSC® TWIN GRASPER®