



A.M.I. HAL-Doppler II System plus A.M.I. RAR Probe and Sleeve suitable for:

> **A.M.I. HAL (Haemorrhoidal Artery Ligation)**

- Minimally-invasive treatment for lower grades of haemorrhoids

> **A.M.I. RAR (Recto Anal Repair)**

RAR is a combination of HAL plus mucopexy

- Treat haemorrhoidal symptoms and prolapsing tissue in one quick session
- No tissue resection necessary
- Low pain and low risk treatment for more severe haemorrhoids

> **Simple, safe and gentle treatment for all grades of haemorrhoids**

A.M.I.® RAR System










Step 1:
HAL -
Haemorrhoidal
Artery Ligation



Step 2a:
RAR -
Recto Anal Repair



Step 2b:
RAR -
Recto Anal Repair

Order Code	A.M.I. RAR System	Technical Details
AHD 204	A.M.I. HAL-Doppler II System Kit consists of: - AHE 203 A.M.I. HAL-II Electronic System incl. HAL-II Adaptor and HAL-II Rechargeable Battery Pack - AHH 001 A.M.I. HAL Handle - RAR 2013 A.M.I. Recto Anal Repair Sleeve - RAR 2012 A.M.I. HAL Probe Single use* - AHK 007 A.M.I. HAL Knotpusher - AHN 006 A.M.I. HAL Needleholder - AHAL 70 A.M.I. HAL Suture - Instructions for use * RAR Probe also available on request	Displays identified artery depth and blood flow velocity. 1 unit with inbuilt printer Unit dimensions 385 x 140 x 260 mm 2 units 2 units 5 units (Dimensions 80 mm, Ø 28 mm) 2 units 2 units 1 box of 36 sutures
AHH 001	A.M.I. HAL Handle  For use with the A.M.I. HAL-II Electronic System to hold the ultrasound probe	1 unit
RAR 2013	A.M.I. Recto Anal Repair Sleeve  Outer sleeve to attach to the A.M.I. HAL Handle	1 unit Delivered non-sterile, steam autoclavable
RAR 2011	A.M.I. Recto Anal Repair Probe  Single use probe for RAR procedure	1 box of 5 probes Delivered sterile, fully disposable
AHS 004	A.M.I. HAL Probe Single use  Single use probe for HAL procedure	1 box of 5 probes Delivered sterile, fully disposable
RAR 2012	A.M.I. Recto Anal Repair Fixation Nut  To attach the A.M.I. RAR Probe or A.M.I. HAL Probe to the A.M.I. HAL Handle	1 unit Delivered non-sterile, steam autoclavable
AHN 006	A.M.I. HAL Needleholder  Specially designed to match the ligation groove inside the HAL Probe	1 unit Delivered non-sterile, steam autoclavable
AHK 007	A.M.I. HAL Knotpusher  To facilitate knot tying inside the A.M.I. Probes	1 unit Delivered non-sterile, steam autoclavable
AHAL 70	A.M.I. HAL Suture  Synthetic, absorbable, 2-0	1 box of 36 sutures 75 cm long, 5/8 circle needle

A.M.I. Headquarters:

A.M.I. GmbH
Im Letten 1
6800 Feldkirch . Austria
t +43 5522 90505-0
f +43 5522 90505-4006
e info@ami.at
www.ami.at